

Molecular Methods in Microbial Ecology

Contact Info: Julie Huber and Nancy Akerman

Schedule: 25 Oct: Introductory Lecture, DNA extraction
Run DNA products on gel


27 Oct: Lecture on PCR
Prepare PCR reactions


1 Nov: Analyze gels from PCR
Lecture on other molecular methods

Readings: Head *et al.* 1998. *Microbial Ecology* 35: 1-21.

Day 1

- Introduction to molecular methods in microbial ecology
- Extract DNA from Winogradsky Columns
- Run DNA on agarose gel






The Challenge for Microbial Ecology

Habitat	Culturability (%)
Seawater	0.001-0.1
Freshwater	0.25
Sediments	0.25
Soil	0.3

**How do you study something you can't
grow in the lab?**

The Solution: Molecular Biology



- Present in all cells- Bacteria, Archaea and Eukaryotes
- Documents of evolutionary history
- Basis of all molecular biological techniques

Environmental Sample

DNA

RNA

What organisms are present?

What genes are present?

What organisms are active?


What genes are being expressed?


Target Genes

Metagenomics


Microarrays

Proteomics





DNA extraction from Winogradsky Columns





DNA Extraction

1. Lyse cell membrane
 - a. Chemically → detergent
 - b. Physically → bead beating
2. Pellet cell membrane, proteins and other cell parts while DNA stays in solution
3. Remove other inhibitors from DNA
4. Mix DNA with acid and salt → stick to filter
5. Wash filter-bound DNA several times with alcohol
6. Elute DNA off membrane with pH 8, low-salt buffer

Choosing a Depth Horizon

- 16S rRNA Bacteria
- 16S rRNA Archaea
- *mcrA* Methanogens
 - Methyl coenzyme M reductase
- *dsrB* Sulfate reducers
 - Dissimilatory bisulfite reductase
- *pmoA* Methanotrophs
 - Particulate methane monooxygenase








Day 1, Part II

- Run an electrophoresis gel of the DNA products extracted from your columns


Basics of Gel Electrophoresis

- The gel is a matrix (like jello with holes)
- DNA is negatively charged- will run to positive
- Smaller fragments run faster than larger ones
- Gel contains Ethidium Bromide, which binds to DNA and fluoresces when hit with UV light
(WEAR GLOVES!!!)




Genomic DNA

The sum total of all DNA from an organism or a community of organisms



What to do

- Mix 10 μl of your DNA with 5 μl loading buffer
- Load in well on gel
- I'll load the ladder
- Run it
- Take a picture of it





Day 2


- Learn about PCR
- Set up PCR reactions using the DNA from your extractions and an assortment of primers

Genomic DNA

The sum total of all DNA from an organism or a community of organisms







Environmental Sample

DNA

RNA

What organisms are present?

What genes are present?

What organisms are active?

What genes are being expressed?

Target Genes

Metagenomics

Microarrays

Proteomics


The Star of the Show: SSU rRNA

- Everybody has it
- Contains both highly conserved and variable regions
 - allows making comparisons between different organisms over long periods of time (evolutionary history)
- Not laterally transferred between organisms
- HUGE and growing database




Ribosomes

- Make proteins
- rRNA is transcribed from rDNA genes




SSU rRNA



Universal Tree of Life

BACTERIA

ARCHAEA



EUKARYA

Modified from Norman Pace

Polymerase Chain Reaction (PCR)

- Rapid, inexpensive and simple way of making millions of copies of a gene starting with very few copies
- Does not require the use of isotopes or toxic chemicals
- It involves preparing the sample DNA and a master mix with primers, followed by detecting reaction products


Every PCR contains:

- A DNA Polymerase (most common, Taq)
- Deoxynucleotide Triphosphates (A, C, T, G)
- Buffer (salt, $MgCl_2$, etc)
- A set of primers, one Forward, one Reverse
- Template DNA


Typical PCR Profile

Temperature	Time	Action
95°C	5 minutes	DNA Taq polymerase activation
35 cycles of: 95°C 54°C 72°C	1 minute 1 minute 1 minute	DNA denaturization Primer annealing Extension creation
72°C	10 minutes	Final extension created

Polymerase Chain Reaction (PCR)




Denature DNA (1 double strand to 2 single strands) and attach primers




Enzyme (DNA Polymerase) creates complementary strands from free nucleotides



Repeat



Geometric increase in PCR product



30 cycles of PCR creates 2^{30} copies of the original template > 1,000,000,000 copies

Things you can optimize

- Temperature and time to activate Taq polymerase
- Temperature and time to allow primer annealing
- Temperature and time for extension
- Concentration of reagents, especially primers, dNTPs, and $MgCl_2$
- Concentration of template DNA
- Number of replication cycles
- Etc...

Beyond 16S

- Identical 16S \neq Identical Function
- Target functional genes


	26J	28U	30U	32U
species 1	TCAAAGATTAAGC-CAT3CATGTCAAAGT--ACAAGCCCCATT A-AG -G7GA-AAACC3C3AATGGCTCAITAAA--TCA			
species 2	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAATCCTCTTGA-GG-GAGA-AACTGC3AAAGGCTCAITAAA--TCA			
species 3	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCCCTCACT A-AG -G7GA-AAACC3C3AATGGCTCAITAAA--TCA			
species 4	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAGGCCG-ATCT- AAG -GCGA-AAACC3C3AATGGCTCAITAAA--TCA			
species 5	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAGGCCG-ATCT- AAG -GCGA-AAACC3C3AATGGCTCAITAAA--TCA			
species 6	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAGGCCG-AACT- AAG -GCGA-AAACC3C3AATGGCTCAITAAA--TCA			
species 7	TGGTGTCTCGTT3CCT3C-TGTCTAAGT--ACAAGCCG-ATTC- AAG -GCGA-AAACC3C3AATGGCTCAITAAA--TCA			
species 8	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCCG-ATTT- AAG -GCGA-AAACC3C3AATGGCTCAITAAA--TCA			
species 9	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCCG-ATGT- AAG -G7GA-AAACC3C3AATGGCTCAITAAA--TCA			
species 10	TCAAAGATTAAGC-CAT3CATGTCCTNNGT--ACA---CCTCTG--GG--GCGA-AAACC3C3AATGGCTCAITAAAA--TCA			
species 11	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCGCTATG--CG--GCGA-AAACC3C3AATGGCTCAITAAA--TCA			
species 12	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCCGCTAGA-CG--GCGA-AAACC3C3AATGGCTCAITAAA--TCA			
species 13	TCAAAGATTAAGC-CAT3CAGGTCCTAAGT--ATAAGCCGAAATA-AA--G7GA-GACC3C3AATGGCTCAITACA--TCA			
species 14	TCAAAGATTAAGC-CAT3CAGGTCCTAAGT--ATGAGCCGAAATA-AAT-G7GA-GACC3C3AATGGCTCAITACA--TCA			
species 15	TCAAAGATTAAGC-CAT3CAGGTCCTAAGT--ACATGCTCTTATA-TATGGTAA-GACT3C3AACCGCTCAITACA--TCG			
species 16	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCCGCTA A-AG -G7GA-AAACC3C3AATGGCTCAITAAA--TCA			
species 17	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCCGCTA A-AG -G7GA-AAACC3C3AATGGCTCAITAAA--TCA			
species 18	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCCGCTA A-AG -G7GA-AAACC3C3AATGGCTCAITAAA--TCA			
species 19	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCCGCTA A-AG -G7GA-AAACC3C3AATGGCTCAITAAA--TCA			
species 20	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCCGCTA A-AG -G7GA-AAACC3C3AATGGCTCAITAAA--TCA			
species 21	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--TCA-ACCTCGTCT-CG--CGGACAACT3C3GATGGCTCAITAAA--TCA			
species 22	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--ACAAGCCGCTA A-AG -G7GA-AAACC3C3AATGGCTCAITAAA--TCA			
species 23	TCAAAGATTAAGC-CAT3CATGTCCTAAGT--TCA TG -G7GA-AAACC3C3AATGGCTCAITAAA--TCA			

16S rDNA

mcrA




A



B

Primers We Are Using

- 16S rRNA Bacteria
- 16S rRNA Archaea
- *mcrA* Methanogens
 - Methyl coenzyme M reductase
- *dsrB* Sulfate reducers
 - Dissimilatory bisulfite reductase
- *pmoA* Methanotrophs
 - Particulate methane monooxygenase



Reagent	Volume (μl) per reaction	# of reactions	final volume
Sterile H ₂ O	26.8		
5X PCR buffer	10		
dNTPs (10mM)	1		
Taq polymerase (5 Units/ μl)	0.2		

Tube	Master mix	Target	Template	Vol	F primer	Vol	R primer	Vol
	μl			μl		μl		μl
1	38	Sulfate reducers	Column DNA	2	dsr1F	5	dsr4R	5
2	38	Methanogens	Column DNA	2	ME1	5	ME2	5
3	38	Methanotrophs	Column DNA	2	pa189F	5	pa682R	5
4	38	Bacteria	Column DNA	2	8F	5	1492R	5
5	38	Archaea	Column DNA	2	21F	5	958R	5
6	38	Archaea	+ control <i>M. jannaschii</i>	2	21F	5	958R	5
7	38	Nothing	- control (water)	2	21F	5	958R	5